



DECISIO

**DECISIO HEALTH PATIENT DASHBOARD
USER GUIDE ADDENDUM – GE MURAL APPLICATION**

DOCUMENT VERSION: LB-20027, REV. B

**DECISIO HEALTH, INC.
520 POST OAK BOULEVARD SUITE 600
HOUSTON, TX 77027**

Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

USING MURAL WITH THE DECISIO PATIENT DASHBOARD

To begin, open the Mural application and login with the user's credentials. If the user encounters the image below, in Figure 1, a user login error has occurred. Please restart the Mural application and try to login again. If the problem persists, contact the System Administrator for resolution.




Figure 1. User Login Error

After successful login, the user will see the patient surveillance screen and the GE logo screen, seen below in Figure 2.



Figure 2. The patient surveillance screen is seen on the left and the GE logo screen on the right.

To view a patient dashboard, the user can choose a patient by clicking on the header of a patient index card in the surveillance screen. The patient dashboard will replace the GE Healthcare logo in the next screen. See figure 2 above with the GE Healthcare logo and Figure 3 below with the GE Healthcare logo replaced by the patient dashboard.



Figure 3. The patient surveillance screen is shown on the left and the patient dashboard is on the right. The patient dashboard has replaced the GE Healthcare logo.

If the user encounters a blurred header and the words “NO MRN” for a patient, the user will be unable to select that patient until the MRN is resolved. Once the MRN is resolved, the blurred header will change back to normal and display the bed number, MRN and patient name. See Figure 4 below. The user can then select the patient to view the dashboard, after the MRN resolved.



Figure 4. Unresolved MRN with header blurred is shown on the left. On the right, once the MRN resolves, the header will return to its normal state.

If the user encounters the image below, in Figure 5, at any time while using the surveillance screen or patient dashboard, a system error has occurred. Please restart the Mural application, login and attempt to access the surveillance screen and patient dashboard again. If the problem persists, contact the System Administrator for resolution.



Figure 5. System error.

Manufactured By:

Decisio Health, Inc.

520 Post Oak Boulevard Suite 600

Houston, TX 77027 www.decisiohealth.com

General Information Email: info@decisiohealth.com

Support Email: support@decisiohealth.com